
|                 |  |
|-----------------|--|
| <b>EN</b><br>5  | <b>Bakri Postpartum Balloon</b><br>Instructions for Use    |
| <b>CS</b><br>7  | <b>Tampónovací balónek Bakri</b><br>Návod k použití        |
| <b>DA</b><br>9  | <b>Bakri postpartum-ballon</b><br>Brugsanvisning           |
| <b>DE</b><br>11 | <b>Bakri Postpartum Ballon</b><br>Gebrauchsanweisung       |
| <b>EL</b><br>14 | <b>Επιλόχειο μπαλόني Bakri</b><br>Οδηγίες χρήσης           |
| <b>ES</b><br>16 | <b>Balón de posparto Bakri</b><br>Instrucciones de uso     |
| <b>FR</b><br>18 | <b>Ballonnet post-partum de Bakri</b><br>Mode d'emploi     |
| <b>HU</b><br>21 | <b>Bakri postpartum ballon</b><br>Használati utasítás      |
| <b>IT</b><br>23 | <b>Palloncino post-parto Bakri</b><br>Istruzioni per l'uso |
| <b>NL</b><br>25 | <b>Bakri-postpartumballon</b><br>Gebruiksaanwijzing        |
| <b>NO</b><br>28 | <b>Bakri-postpartumballong</b><br>Bruksanvisning           |
| <b>PL</b><br>30 | <b>Balon poporodowy Bakriego</b><br>Instrukcja użycia      |
| <b>PT</b><br>32 | <b>Balão pós-parto Bakri</b><br>Instruções de utilização   |
| <b>SV</b><br>35 | <b>Bakri postpartumballong</b><br>Bruksanvisning           |

### Improper Placement

Nesprávné umístění  
Ukorrekt anläggelse  
Falsche Platzierung  
Εσφαλμένη τοποθέτηση  
Colocación incorrecta  
Mise en place incorrecte  
Nem megfelelő behelyezés  
Posizionamento errato  
Verkeerd geplaatst  
Feil plassering  
Umieszczenie nieprawidłowe  
Colocação incorreta  
Felaktig placering




Figure 1

### Proper Placement


Správné umístění  
Korrekt anläggelse  
Richtige Platzierung  
Σωστή τοποθέτηση  
Colocación correcta  
Mise en place correcte  
Megfelelő behelyezés  
Posizionamento corretto  
Correct geplaatst  
Riktig plassering  
Umieszczenie prawidłowe  
Colocação correta  
Korrekt placering


**Figure 2**


**Figure 3**


**Figure 4**


**Figure 5**

**Figure 6**

**Figure 7**

**Figure 8**

## BAKRI POSTPARTUM BALLOON

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### INTENDED USE

Use of this device is intended to provide temporary control or reduction of postpartum uterine bleeding when conservative management is warranted.

### CONTRAINDICATIONS

- Arterial bleeding requiring surgical exploration or angiographic embolization
- Cases indicating hysterectomy
- Pregnancy
- Cervical cancer
- Purulent infections in the vagina, cervix, or uterus
- Untreated uterine anomaly
- Disseminated intravascular coagulation
- A surgical site which would prohibit the device from effectively controlling bleeding

### WARNINGS

- This device is intended as a temporary means of establishing hemostasis in cases indicating conservative management of postpartum uterine bleeding.
- The device should not be left indwelling for more than 24 hours.
- Always inflate the balloon with sterile liquid. Never inflate with air, carbon dioxide or any other gas.
- The maximum inflation is 500 ml. Do not overinflate the balloon.
- Clinical data to support the safety and effectiveness of the Bakri Postpartum Balloon in the setting of uterine atony are limited. Patients in whom this device is being used should be closely monitored for signs of worsening bleeding and/or disseminated intravascular coagulation (DIC). In such cases, emergency intervention per hospital protocol should be followed.
- There are no clinical data to support use of this device in the setting of DIC.
- Patient monitoring is an integral part of managing postpartum hemorrhage. Signs of deteriorating or non-improving condition should lead to a more aggressive treatment and management of patient uterine bleeding.
- Patient urine output should be monitored while the Bakri Postpartum Balloon is in use.

### PRECAUTIONS

- Avoid excessive force when inserting the balloon into the uterus.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

**IMPORTANT:** In order to maximize tamponade effect, constant pressure must be effected between the balloon wall and the tissue surface. Pressure can be achieved and maintained by applying gentle traction to the balloon shaft, before securing to the patient's leg, or attaching to a weight, not to exceed 500 grams.

### INSTRUCTIONS FOR USE

#### Transvaginal Placement

1. Determine uterine volume by direct examination or ultrasound examination.

2. Insert the balloon portion of the catheter into the uterus, making certain that the entire balloon is inserted past the cervical canal and internal ostium.
3. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.

### **Transabdominal Placement, Post Cesarean Section**

1. Determine uterine volume by direct examination (intraoperative) or ultrasound examination (postoperative).
2. From above, via access of the Cesarean incision, pass the tamponade balloon, inflation port first, through the uterus and cervix. **NOTE:** Stopcock may be removed to aid placement and reattached prior to filling balloon.
3. Have an assistant pull the shaft of the balloon through the vaginal canal until the deflated balloon base comes into contact with the internal cervical ostium.
4. Close the incision per normal procedure, taking care to avoid puncturing the balloon while suturing.  
**NOTE:** Ensure that all product components are intact prior to inflating the balloon.

### **Balloon Inflation with Syringe**

**WARNING:** Always inflate the balloon with sterile liquid. Never inflate with air, carbon dioxide or any other gas.

**WARNING:** The maximum inflation is 500 ml. Do not overinflate the balloon.

**NOTE:** To ensure that the balloon is filled to the desired volume, it is recommended that the predetermined volume of fluid be placed in a separate container, rather than relying on a syringe count to verify the amount of fluid that has been instilled into the balloon.

1. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.
2. Using the enclosed syringe, begin filling the balloon to the predetermined volume through the stopcock.
3. Once the balloon has been inflated to the predetermined volume, confirm placement via ultrasound.  
**(Fig. 1) NOTE:** See **Fig. 1** for proper placement.
4. Apply gentle traction to the balloon shaft to ensure proper contact between the balloon and tissue surface. To maintain tension, secure the balloon shaft to the patient's leg or attach to a weight, not to exceed 500 grams.  
**NOTE:** To maximize the tamponade effect, counter pressure can be applied by packing the vaginal canal with iodine- or antibiotic-soaked vaginal gauze.
5. Connect the drainage port to a fluid collection bag to monitor hemostasis. **NOTE:** To adequately monitor hemostasis, the balloon drainage port and tubing should be flushed clear of clots with sterile isotonic saline.
6. Monitor the patient continuously for signs of increased bleeding and uterine cramping.

### **Balloon Inflation with Rapid Instillation Components**

See pictorial instructions, **Figs. 2-8**, at the front of this booklet.

**NOTE:** Ultrasound should be used to confirm proper placement of the balloon once the balloon is inflated to the predetermined volume.

### **Balloon Removal**

**NOTE:** The maximum indwell time is 24 hours. The balloon may be removed sooner upon the physician's determination of hemostasis.

1. Remove tension from the balloon shaft.
2. Remove any vaginal packing.
3. Using an appropriate syringe, aspirate the contents of the balloon until fully deflated.
4. Gently retract the balloon from the uterus and vaginal canal and discard.
5. Monitor patient for signs of bleeding.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

## REFERENCE

Y. N. Bakri, A. Amri, F. Abdul Jabbar: "Tamponade balloon for obstetrical bleeding," *International Journal of Gynecology and Obstetrics*, 74 (2001) 139-142.

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© COPYRIGHT COOK 2012
**COOK IRELAND LTD.**